

# KlenTaq DNA Polymerase

Cat # : PT-GL-KTAQ      Size : 250U (5U/ $\mu$ L)  
         SA-GL-KTAQ      Size : 25U (5U/ $\mu$ L)

## Product Description:

KlenTaq DNA Polymerase is missing the N-terminal portion of the gene, encoding *Thermus aquaticus* (*Taq*) DNA polymerase, leaving a highly active and even more thermal stable DNA polymerase activity. This enzyme keeps significant activity after exposure to 99°C or repeated exposure to 98°C. The missing N-terminal portion is homologous to the 5'→3' exonuclease region of *E. coli*. DNA polymerase I. KlenTaq is therefore similar to, yet distinct from, Hoffman LaRoche's Stoffel Fragment.

## Components: Store at -20 °C (Store buffer at 4 °C)

- \* 1 tube of 50 $\mu$ L KlenTaq DNA polymerase in KlenTaq Storage Buffer (SA: 5uL)
- \* 1 tube 10X KlenTaq PCR Buffer 1mL (SA: 100uL)

## KlenTaq Storage Buffer :

50% glycerol (v/v), 50 mM (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>, 20 mM Tris-HCl pH8.55, 0.1 mM EDTA, 10 mM mercaptoethanol, no gelatin, and 0.5% Thesit.

- \* Thesit (from Boehringer-Mannheim) with no OD<sub>280</sub> is a superior nonionic detergent for *Taq* polymerase storage. It replaces the previously used triton X100, NP40, and/or Tween-20.

## 10X KlenTaq PCR Buffer :

50 mM Tris-HCl pH9.1, 1.6 mM (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>, 3.5 mM MgCl<sub>2</sub>, and 150  $\mu$ g/mL BSA.

- \* The dNTP concentration can vary from 50  $\mu$ M each to 1.2 mM each, but 200  $\mu$ M is normally used.
- \* Suggested enzyme amount is 0.5  $\mu$ L per 100  $\mu$ L reaction.

The amount of KlenTaq needed to perform a reaction depends on the length of incorporation. For 100  $\mu$ L reaction, 1  $\mu$ L KlenTaq will catalyze 15-20 reactions of 500 bp and 8-10 reactions of 1 kb and 3-5 reactions of 2 kb template DNA. Since excess KlenTaq is harmless, we conservatively recommend 0.5  $\mu$ L per reaction in order to be sure everything up to 2.5 kb will work.

- \* Please see **General Protocol** in the next page.

## General Protocol:

Reaction Mix Preparation Table

Component	Volume
10X KlenTaq PCR Buffer	5 $\mu$ L
dNTP (2.5mM each)	4 $\mu$ L (200 $\mu$ M)
Upstream primer (10 $\mu$ M)	1 $\mu$ L
Downstream primer (10 $\mu$ M)	1 $\mu$ L
DNA template	X* $\mu$ L
KlenTaq DNA polymerase	0.5 $\mu$ L
Sterile H <sub>2</sub> O	Up to 50 $\mu$ L

\* X is variable depend on user's condition.

Suggested PCR Program

Stage	Temperature	Time
Initial denature	94-95°C	1-5 min
Denature	25-40 cycles	94-95°C
Annealing		50-60°C*
Extension		72°C
Final extension	72°C	5-10 min
Preservation	4-10°C	

\* Annealing temperature is depending on the T<sub>m</sub> value of primers.

\* KlenTaq DNA Polymerase is more thermal stable than normal *Taq* DNA Polymerase. Therefore, denature temperature can be slightly arisen according to the need of experiment.

- For experimental sample, please test for your own cycle conditions.
- This enzyme is stored in -20 °C for at least 1 year.

\*\*\*Research Use Only\*\*\*

Please do not hesitate to contact us while you have any questions.  
Manufactured for and distributed by Protech Technology Enterprise Co., Ltd.  
Tel : +886-2-2655-7677 / FAX : +886-2-2655-7601 / Toll Free : 0800-231-530  
E-MAIL: [service@bio-protech.com.tw](mailto:service@bio-protech.com.tw) ; [tech@bio-protech.com.tw](mailto:tech@bio-protech.com.tw)  
Web Site: <http://www.bio-protech.com.tw>