

Date: 2011.05

Version 11_01

4X Upper Gel Buffer

Cat No : PT-4X-UP-1L Size : 1 L

Description :

4X Upper Gel Buffer contain Tris and SDS at the appropriate *pH* and concentration for convenient and quick preparation of polyacrylamide gels for protein analysis by SDS-PAGE. Dilution of the Stacking Gel 4X buffer (0.5M Tris *pH* 6.8 and 0.4% SDS) to a 1X solution results in the standard, 0.125M Tris, 0.1% SDS concentrations routinely used in the preparation of the stacking gel for SDS-PAGE.

Component : *Store at RT*

4X Upper Gel Buffer (0.5M Tris *pH* 6.8 and 0.4% SDS)

Research Use Only

Please do not hesitate to contact us while you have any questions.

Manufactured for and distributed by Protech Technology Enterprise Co., Ltd.

TEL : +886-2-2655-7677 / FAX : +886-2-2655-7601 / Toll Free : 0800-231-530

E-MAIL: service@bio-protech.com.tw ; tech@bio-protech.com.tw

WEB: <http://www.bio-protech.com.tw>