

Pro Taq Plus DNA Polymerase

Certificate of Analysis

Cat No : PT-526-2U

Size : 500 units

Concentration : 2 units/ μ l

Volume : 250 μ l

Unit Definition :

One unit is the amount of enzyme that will convert 10 nmoles of dNTPs into acid-insoluble form in 30 minutes at 70°C under that stated assay.

Assay Conditions :

10 mM Tris-HCl (pH 9.0)(25°C), 50 mM KCl, 1.5 mM MgCl₂, 0.1% (w/v) gelatin, 1% Triton X-100 incubate at 72°C

Storage Buffer :

50 mM Tris-HCl (pH 8.0), 1 mM EDTA, 1 mM DTT, 50% (v/v) glycerol.

Reaction Buffer :

10X Buffer contains : 100 mM Tris-HCl (pH 9.0), 500 mM KCl, 0.1% (w/v) gelatin, 15 mM MgCl₂, 1% Triton X-100

Endonuclease Assay :

Incubation of 10 units of the enzyme in assay buffer with 0.5 μ g of pBR322 DNA for one hour at 70°C give < 5% conversion of RFI to RF II DNA.

Exonuclease Assay :

Incubation of 20 units of the enzyme with 1 μ g lambda DNA for 16 hours at 65°C in the *Tth* I restriction enzyme buffer does not produce any detectable degradation of the DNA.

For Research Using Only.

Please do not hesitate to contact us if you have any questions.

Manufactured for and distributed by Protech Technology Enterprise Co.,Ltd

TEL : +886-2-2655-7677 / FAX : +886-2-2655-7601 / Toll Free : 0800-231-530

E-MAIL: service@bio-protech.com.tw ; tech@bio-protech.com.tw

WEB: <http://www.bio-protech.com.tw>