

M13 Forward Primer (-20)

Certificate of Analysis

Catalog number : N520-02

Quantity : 2 µg

Volume : 200 µl

Concentration : 10 ng/µl (1.62 µM)

Description :

Universal primer for sequencing from most M13 forward primer sites (-20), 16 mers. M13 forward primer is purified by RP-HPLC.

Sequence :

5'-g TAA AAC gAC ggC CAg-3'

Storage condition :

The solutions are stable at -20°C. However, for long term storage, we recommend storing in aliquots at -70°C.

For Research Using Only.

Please do not hesitate to contact us if you have any questions.

Manufactured for and distributed by Protech Technology Enterprise Co.,Ltd

TEL : +886-2-2655-7677 / FAX : +886-2-2655-7601 / Toll Free : 0800-231-530

E-MAIL: service@bio-protech.com.tw ; tech@bio-protech.com.tw

WEB: <http://www.bio-protech.com.tw>