

MSPS SPECIFICATION

Model	<i>MSPS-200</i>	<i>MSPS-10</i>
MicroTip	200ul reference tip	10ul reference tip
Max. Volume	1500 tips	3000 tips
Capacity (for Protech Racks)	(8 x 96well)racks / bench	(8 x 96well)racks / bench
	(8 x 96well)plates / bench	(8 x 96well)plates / bench
Operation Time	100 sec per rack	100 sec per rack
Operation System	FBS-20MCT Control system	FBS-20MCT Control system
Power Supply	AC220V, 50/60Hz, 100W	AC220V, 50/60Hz, 100W
Size (L×W×H)	95cm × 85cm × 180cm	95cm × 85cm × 180cm
Net Weight	300kg	300kg
Abstract	1 x Base 1 x XY Platform 1 x Carrying Platform 1 x Vibration Tray 1 x Auto Discharge Mechanism 1 x Control Unit	
Brand	Protech (Made in Taiwan)	