

# Oligo dT<sub>12-18</sub> Primer

## *Certificate of Analysis*

**Catalog number :** N420-01

**Quantity :** 25 µg

**Concentration :** 1 µg/µl (≈ 0.2 mM)

**Volume :** 25 µl

**Description :**

The oligo dT primer consists 12-18 dT bases and is designed to primer poly A+ RNA for first-strand cDNA synthesis. Oligo dT<sub>12-18</sub> primer is purified by RP-HPLC.

**Storage condition :**

The solutions are stable at -20°C. However, for long term storage, we recommend storing in aliquots at -70°C.

*For Research Using Only.*

*Please do not hesitate to contact us if you have any questions.*

*Manufactured for and distributed by Protech Technology Enterprise Co.,Ltd*

TEL : +886-2-2655-7677 / FAX : +886-2-2655-7601 / Toll Free : 0800-231-530

E-MAIL: service@bio-protech.com.tw ; tech@bio-protech.com.tw

WEB: <http://www.bio-protech.com.tw>