

Random Primer

Certificate of Analysis

Catalog number : N410-01

Quantity : 25 µg

Concentration : 1 µg/µl (≈ 0.5 mM)

Description :

The random primer consists 6 N bases and is designed to primer RNA for first-strand cDNA synthesis. Random primer purified by RP-HPLC.

Storage condition :

The solutions are stable at -20°C. However, for long term storage, we recommend storing in aliquots at -70°C.

For Research Using Only.

Please do not hesitate to contact us if you have any questions.

Manufactured for and distributed by Protech Technology Enterprise Co.,Ltd

TEL : +886-2-2655-7677 / FAX : +886-2-2655-7601 / Toll Free : 0800-231-530

E-MAIL: service@bio-protech.com.tw ; tech@bio-protech.com.tw

WEB: <http://www.bio-protech.com.tw>