

SP6 Promoter Primer

Certificate of Analysis

Catalog number : 69349-1

Quantity : 2 μg

Concentration : 10 $\text{ng}/\mu\text{l}$ (1.62 μM)

Volume : 200 μl

Description :

Universal primer for sequencing from most SP6 primer sites, 19 mers. SP6 promoter primer is purified by RP-HPLC.

Sequence :

5'-gAT TTA ggT gAC ACT ATA g-3'

Storage condition :

The solutions are stable at -20°C . However, for long term storage, we recommend storing in aliquots at -70°C .

For Research Using Only.

Please do not hesitate to contact us if you have any questions.

Manufactured for and distributed by Protech Technology Enterprise Co.,Ltd

TEL : +886-2-2655-7677 / FAX : +886-2-2655-7601 / Toll Free : 0800-231-530

E-MAIL: service@bio-protech.com.tw ; tech@bio-protech.com.tw

WEB: <http://www.bio-protech.com.tw>