

# T7 Promoter Primer

## *Certificate of Analysis*

**Catalog number :** 69348-1

**Quantity :** 2 µg

**Concentration :** 10 ng/µl (1.64 µM)

**Volume :** 200 µl

**Description :**

Universal primer for sequencing from most T7 primer sites, 20 mers. T7 promoter primer purified by RP-HPLC.

**Sequence :**

5'-TAA TAC gAC TCA CTA TAg gg-3'

**Storage condition :**

The solutions are stable at -20°C. However, for long term storage, we recommend storing in aliquots at -70°C.

*For Research Using Only.*

*Please do not hesitate to contact us if you have any questions.*

*Manufactured for and distributed by Protech Technology Enterprise Co.,Ltd*

TEL : +886-2-2655-7677 / FAX : +886-2-2655-7601 / Toll Free : 0800-231-530

E-MAIL: service@bio-protech.com.tw ; tech@bio-protech.com.tw

WEB: <http://www.bio-protech.com.tw>