

# Deoxynucleotide Triphosphates (dNTP) Pre-Mixed

## *Certificate of Analysis*

**Catalog number :** PTM8010

**Quantity :** 1 ml

**Concentration :** 10 mM for each of dATP, dCTP, dGTP and dTTP

**Description :**

The dNTP is made by dissolving ultra-pure triphosphate powders in buffer at pH 7.4. The concentration was quantitated by spectral analysis and then adjusted to a final concentration of 10 mM

**Storage condition :**

The solutions are stable at -20°C. However, for long term storage, we recommend storing in aliquots at -70°C.

***For Research Using Only.***

***Please do not hesitate to contact us if you have any questions.***

*Manufactured for and distributed by Protech Technology Enterprise Co.,Ltd*

TEL : +886-2-2655-7677 / FAX : +886-2-2655-7601 / Toll Free : 0800-231-530

E-MAIL: service@bio-protech.com.tw ; tech@bio-protech.com.tw

WEB: <http://www.bio-protech.com.tw>