

Deoxynucleotide Triphosphates (dNTP) Pre-Mixed

Certificate of Analysis

Catalog number : PTM8014

Quantity : 1 ml

Concentration : 2.5 mM for each of dATP, dCTP, dGTP and dTTP

Description :

The dNTP is made by dissolving ultra-pure triphosphate powders in buffer at pH 7.4. The concentration was quantitated by spectral analysis and then adjusted to a final concentration of 2.5 mM

Storage condition :

The solutions are stable at -20°C. However, for long term storage, we recommend storing in aliquots at -70°C.

For Research Using Only.

Please do not hesitate to contact us if you have any questions.

Manufactured for and distributed by Protech Technology Enterprise Co.,Ltd

TEL : +886-2-2655-7677 / FAX : +886-2-2655-7601 / Toll Free : 0800-231-530

E-MAIL: service@bio-protech.com.tw ; tech@bio-protech.com.tw

WEB: <http://www.bio-protech.com.tw>